

# Ceglu™ coating solution

500 mL x 1 bottles / box

Ceglu™ is a chemically defined, animal-origin-free scaffold material with excellent homogeneity and storage stability. Ceglu™ coating solution can provide a uniform surface on variety of culture substrates such as plates, flasks and non-woven fabrics to support optimal stem cell culture.

## Properties

- Storage: Store at 1-30 °C.
- Shelf Life: Stable until expiry date on label.

## Precautions

- See SDS:  
**NO OPEN FLAMES ALLOWED. This product contains hazardous liquids.**  
**Please handle with care, using protective equipment in appropriate conditions.**  
**Please refer to the Safety Data Sheet (SDS) before use.**



- Single Use Only:  
This product is designed for single use only. Do not reuse.
- Intended Use:  
This product is for research purposes only and should not be used for therapeutic or diagnostic purposes. The product may not be sold to a third party or used for commercial purposes without our permission.

## Key Features

- Safety:  
Ceglu™ is a fully chemically synthesized scaffold and contains no animal-derived materials. It is also a hydrophobic polymer, which reduces the risk of leaching during culture. This product is Non-pyrogenic.
- Stability:  
It can be stored for long periods at room temperature.
- Culture Quality:  
This product enables attachment and proliferation of stem cells and maintenance of undifferentiation.

## How to use

- Contact us at [support\\_life@sekisui.com](mailto:support_life@sekisui.com)



For more information about Ceglu™ and our other products, please visit our website.



**SEKISUI CHEMICAL CO., LTD.**

2-10-4 Toranomon, Minato-ku, Tokyo, Japan  
Contact : [support\\_life@sekisui.com](mailto:support_life@sekisui.com)